

## **Appointment 1 Consultation and Diagnostic Models**

Patient education and consultation is the first step in any prosthetic procedure. After treatment has been discussed, full arch diagnostic models and accurate bite registrations are taken for laboratory evaluation. Send your case to ROE requesting custom trays. Your case will be carefully evaluated and the multitude of alternative attachments will be considered. We will return custom trays for the final impression. We will also include the suggested treatment plan and estimate of laboratory fee.

## **Appointment 2 Final Impressions**

It is imperative that you select an impression material that can be accurately poured a second time. Utilizing the laboratory constructed custom trays take a full arch upper and lower impression. Master impressions must include all landmarks necessary for partial denture construction. Send master impression (unpoured), opposing model, study models, and detailed instructions to ROE.

We will perform the preparatory work constructing one model with removable dies for the crown and bridge construction and a second full arch model with dies poured in a resin which will be used for the removable appliance. If you prefer to fabricate your own crown and bridge model, please send the master impression to the laboratory for the pouring of the partial model, stabilized bite blocks will be made and returned in the absence of definitive centric stop. If definitive centric stops are present the case can be mounted and case construction can begin, eliminating appointment three.

## **Appointment 3 Registrations**

Centric and esthetic registrations are taken with the stabilized biteblocks provided. Please include as much cosmetic information as possible. At this phase of construction the majority of the laboratory work will be completed. Crowns will be constructed, the partial will be cast, setup and returned for try-in. Laboratory time requirements between appointment 3 and 4 can be up to three weeks depending upon precision attachment availability.

## **Appointment 4 Try-in**

Try-in the case to verify fit, function and esthetics. Make any necessary adjustments and return the case to ROE for completion.

## **Appointment 5 Insertion**

The completed case is inserted. Centric, protrusive and lateral excursions are checked and another ecstatic patient is dismissed. We suggest the patient return to the office for a post insertion bite check in 48 to 72 hours.