

6 Simple Steps to a Perfectly Placed Implant




1

Radiographic template is produced by the laboratory



2

CT Scanning of patient and template



3

CT scan diagnosis and implant planning using Scan2Guide software



4

S2G software creates virtual surgical template



5

Utilizing rapid manufacturing ROE produces surgical iGuide template




6

Simplified and precise placement of implant using iGuide template

Scan2Guide™

- Simple and precise planning of implant surgeries



iGuide™

- Automated process
- Accurate implant placement
- Flapless surgery
- Decreased healing times
- Predictable prosthetic outcomes

iDent™
Simple Yet Sophisticated

9565 Midwest Ave. Garfield Hts, OH 44125
800 228 6663 ext 314
alan@roedentallab.com
www.roedentallab.com