

*Atlantis manufactures abutments in Titanium or Atlantis GoldHue<sup>®</sup> compatible with the implant systems and sizes listed below. Please see other side for compatibility in the Atlantis<sup>®</sup> Abutment in Zirconia<sup>™</sup>. Please consult the Atlantis website or contact Atlantis Customer Service for the most current implant compatibility.*

Implant Manufacturer	Implant Type & Diameter	Abutment Screw Interface Type
<b>Internal Connection Interfaces</b>		
Astra Tech	OsseoSpeed <sup>™</sup> 4.0 <sup>1</sup> , 3.5/4.0, 4.5/ 5.0	0.050" Hex Drive
Biomet 3i	Certain <sup>®</sup> Microminiplant <sup>™</sup> 3.4/3.25 Certain <sup>®</sup> 3.75/4.0, 5.0, 6.0, 4/5, 5/6	0.048" Hex Drive
Nobel Biocare <sup>®</sup>	Replace <sup>®</sup> NP(Narrow Platform 3.5) Replace RP (Regular Platform 4.3) Replace WP (Wide Platform 5.0, 6.0)	Nobel Biocare Unigrip <sup>®</sup>
Straumann <sup>™</sup>	Standard/Standard Plus Implant RN (Regular Neck) 4.8	Straumann Drive Compatible
Zimmer Dental <sup>®</sup>	Tapered Screw-Vent <sup>®</sup> 3.5, 4.5, 5.7 Screw-Vent <sup>®</sup> 3.5, 4.5 SwissPlus <sup>®</sup> 4.8	0.050" Hex Drive or 1.25mm Hex Drive
<b>External Hex Interfaces</b>		
BioHorizons <sup>®</sup>	3.5, 4.0, 5.0	0.050" Hex Drive
Biomet 3i <sup>™</sup>	Microminiplant <sup>™</sup> 3.4/3.25 Miniplant <sup>®</sup> 4.1/3.25 Standard 3.75, 4.0 Wide 5.0, 6.0 XP 3/4 , 4/5 , 5/6	0.048" Hex Drive
Innova Endopore <sup>®</sup> Innova Entegra <sup>™</sup>	3.5, 4.1, 5.0	0.050" Hex Drive
Lifecore <sup>®</sup>	Small Diameter 3.3, 3.4 Regular Diameter 3.75, 4.0 Wide Diameter 5.0, 5.5, 6.0	0.050" Hex Drive
Nobel Biocare <sup>®</sup>	Branemark <sup>®</sup> NP (Narrow Platform 3.3) Branemark <sup>™</sup> RP (Regular Platform 3.75, 4.0) Branemark <sup>™</sup> WP (Wide Platform 5.0, 5.5) SteriOss Replace <sup>®</sup> 3.1, 3.5, 4.3, 5.0, 6.0 SteriOss HL <sup>™</sup> 3.25, 3.8, 4.5, 5.0, 6.0	Nobel Biocare Unigrip <sup>®</sup>
Sterngold-ImplaMed <sup>®</sup>	3.3, 3.75, 4.0, 5.0	0.050" Hex Drive
Zimmer Dental <sup>®</sup>	Spline <sup>®</sup> 3.25, 3.75/4.0, 5.0	0.050" Hex Drive or 1.25mm Hex Drive

All orders for Atlantis Abutments are shipped with an Atlantis HT<sup>™</sup> Titanium Abutment Screw at no extra charge. These abutment screws meet or exceed the implant manufacturer's recommendation for torque. A retaining screw manufactured by Zimmer Dental is included with friction-fit internal hex and Spline<sup>®</sup> implant abutment orders at no extra charge. The implant manufacturer's recommended torque should always be used when placing the abutment screw. Please consult the Atlantis *Instructions for Use* for more information.

<sup>1</sup> Astra Tech "previous" 4.0 with packaging color-coded blue.



Atlantis manufactures Atlantis Abutments in Zirconia™ compatible with the implant systems and sizes listed below. Please see other side for compatibility in Atlantis Abutments in Titanium and Atlantis GoldHue® Abutments. Please consult the Atlantis website or contact Atlantis Customer Service for the most current implant compatibility.

Implant Manufacturer	Implant Type & Diameter	Abutment Screw Interface Type
<b>Internal Connection Interfaces</b>		
Astra Tech	OsseoSpeed™ 3.5/4.0, 4.5/ 5.0	0.050" Hex Drive
Biomet 3i	Certain® Microminiplant™ 3.4/3.25 Certain® 3.75, 4.0, 5.0, 6.0, 4/5, 5/6	0.048" Hex Drive
Nobel Biocare	Replace® RP (Regular Platform 4.3) Replace WP (Wide Platform 5.0, 6.0)	Nobel Biocare Unigrip®
Zimmer Dental®	Tapered Screw-Vent® 3.5, 4.5, 5.7 Screw-Vent® 3.5, 4.5	0.050" Hex Drive or 1.25mm Hex Drive
<b>External Hex Interfaces</b>		
BioHorizons®	4.0, 5.0	0.050" Hex Drive
Biomet 3i™	Miniplant® 4.1/3.25 Standard 3.75, 4.0 Wide 5.0, 6.0 XP 3/4 , 4/5 , 5/6	0.048" Hex Drive
Innova Endopore® Innova Entegra™	4.1, 5.0	0.050" Hex Drive
Lifecore®	Regular Diameter 3.75, 4.0 Wide Diameter 5.0, 5.5, 6.0	0.050" Hex Drive
Nobel Biocare®	Branemark™ RP (Regular Platform 3.75, 4.0) Branemark™ WP (Wide Platform 5.0, 5.5) SteriOss Replace® 4.3, 5.0, 6.0 SteriOss HL™ 3.8, 4.5, 5.0, 6.0	Nobel Biocare Unigrip®
Sterngold-ImplaMed®	3.75, 4.0, 5.0	0.050" Hex Drive

All orders for Atlantis Abutments in Zirconia are shipped with an Atlantis HTZ™ Titanium Abutment Screw at no extra charge. These abutment screws, which have been specifically designed for use with Atlantis Abutments in Zirconia, meet or exceed the implant manufacturer's recommendation for torque. *The implant manufacturer's recommended torque should always be used when placing the abutment screw.* Please consult the Atlantis Instructions for Use for more information.



877-828-5268  
Call Atlantis

support@atlantiscomp.com  
Email Atlantis

www.atlantiscomp.com  
Learn About Atlantis