



9565 Midwest Ave. Garfield Hts., Ohio 44125
 Local 216 663 2233 Toll Free 800 228 6663
 Fax 216 663 2237 Martins Ferry 740 633 9401
 www.roedentallab.com

Dr. Name _____ Phone _____

Patient Name _____ Age _____ Sex M F

Return By 5:00 P.M. On _____ Send Rx Shipping Labels Boxes
 Please Request Case One Day Prior To Appointment Date

FIXED RESTORATIONS

RESTORATION/TOOTH #

- Crown _____
- Bridge _____
- Pontic _____
- Inlay _____
- Veneer _____
- Post & Core _____

MATERIAL

- Porcelain to Metal
- Premium Porc. to Metal
- Procera
- Pressable Ceramic
- Belleglass w/Vectris
- Full Cast
- Acrylic Provisionals

RETURN

- Evaluate & Call
- Dr. Trim Dies
- Diagnostic Waxup
- Cast Try-in
- Bisque Bake
- ★Complete

ALLOY – PORCELAIN TO METAL

- Captek – Yellow (88%Au)
- Hi Noble – Yellow (73%Au) (Premium Porc.)
- Hi Noble – White (67%Au)
- Non-Precious – White
- Titanium – White (Full Cov./Single Crn.)

METAL DESIGN

- ★Lingual Band _____
- Porc. Butt Margin _____
- 360° Porc. On Metal _____
- 360° Metal Band _____
 - 1mm Buccal
 - 1/2mm Buccal
- Metal Occlusal _____
- Metal Lingual _____

TOOTH #

ALLOY – FULL CAST

- Hi Noble – Yellow (77%Au)
- ★Hi Noble – Yellow (56%Au)
- Hi Noble – White (67%Au)

OTHER SERVICES

- Soft Tissue Model
- Crown Under Partial
- Crown Designed For Partial
- Twin Clasp Technique

PONTIC DESIGN




Implant and precision attachment services are available. We encourage study model and treatment plan be submitted for design and estimate of laboratory fees.

Lab Use Only

Dr. called by _____ Date _____

SHADE MAP / CHARACTERIZATION / CONTOUR

Shade _____ Stump Shade _____ (Required for all ceramics)



OCCLUSAL STAIN ★Light Medium Dark None

TRANSLUCENCY ★Standard Increased Minimal Drawn on Shade Map Above

STUDY MODEL FOR Guidance Exact Duplication Incisal Edge Position

CONTOUR ★Match Adjacent/Partner Match Study Model Ideal

OCCLUSION ★Light Full Occlusion Out of Occlusion

IF DEFICIENT OCCLUSAL SPACE Reduction Coping Reduce/MarkOpp. Call


CONTACTS ★Standard Broad Point

REASON FOR RESTORATION Close Spaces Discoloration Aesthetics

FUTURE RESTORATIONS PLANNED _____

Instructions

If no selection is indicated ROE Standards (★) will be used.



Signature _____ License # _____

The person signing this work order accepts responsibility for payment and agrees to pay all collection costs including attorney's fees. A 1 1/2% (18%/yr.) finance charge will be added to all balances due over 30 days.