

# ROE Dental Lab Densply to Ivoclar Mould Conversion

If no mould is provided the following will be selected:

Male	Female
Small Arch = A22	Small Arch = A21
Medium Arch = A25	Medium Arch = A42
Large Arch = A17	Large Arch = A37

If a Dentsply mould is provided it will be converted as follows:

<b>TAPERING</b>	
41J	A17
42D	A41
42F	A22
42G	A66
43D	A11
43F	A13
45H	A26
<b>TAPERING OVOID</b>	
52C	A21
54F	A26
55D	A13
55F	A24B
<b>OVOID</b>	
62D	A22
62G	A37
62G	A37
<b>SQUARE TAPERING OVIID</b>	
74E	A24
75E	A24B
75G	A37
A84	A24
<b>SQUARE</b>	
11G	A25
11H	A17
12E	A66
12G	A15 cent, A26 lat, A15 cusp
13E	A41
<b>Square Tapering</b>	
21C	A42 cent, A41 lat, A42 cusp
21D	A41
21J	A17
21X	A25
22C	A22 cent, A21 lat, A21 cusp
22E	A24B
22G	A34
22H	Dentsply teeth should be used
24F	Dentsply teeth should be used
25G	A34
<b>Square Ovoid</b>	
31F	A25 cent, A25 lat, A24B cusp
32B	A21
32E	A41
32G	A37