

Lab Preferences of Dr. John Doe

Doctor Name

John Doe

Email

roetestemail@protonmail.com

Case Number

123456

General Fixed Preferences

1. Occlusal Contacts

*Light occlusion, slightly out of occlusion (ROE standard)

2. Interproximal Contacts

*Standard contacts - contact with the adjacent teeth (ROE standard)

3. Interproximal Contact Shape

*Standard shape (ROE standard)

4. Model free?

*Yes, whenever possible (ROE standard)

5. Occlusion with a close bite

*Trim opposing up to 0.75mm, mark in black, and send note (ROE standard)

6. Post & Core Copings

*Post and core as unit (ROE standard)

7. Bridge Pontic Design

*Modified ridgelap (ROE standard)

8. Unclear Margins

*Call office (ROE standard)

9. Adjacent Tooth Undercut

*Call office (ROE standard)

10. No Bite Enclosed

*Hand mount if possible (ROE standard)

11. Die Has Undercuts

*Block out undercuts (ROE standards)

Fixed Implant

12. Anterior Lingual Margins

*Crest of tissue (ROE standard)

13. Labial Margins

*1.5mm subgingival (ROE standard)

14. Posterior Lingual Margins

*Crest of tissue (ROE standard)

15. Buccal Margins

*1.0mm subgingival (ROE standard)

16. Interproximal Margins

*0.75mm subgingival (ROE standard)

17. Soft Tissue

*Standard tissue blanching - up to 1 mm of tissue adjustment (ROE standard)

18. Implant Abutment & Restoration Design

*Custom Milled Abutment - Cement Retained (no hole in crown) (ROE standard)

19. Always add pink to the Printed Try-In's?

*No, see Rx (ROE standard)

20. Email full arch case design review?

*No, follow my instructions (ROE standard)

21. Parts

*ROE non-OEM parts preferred but use OEM parts when necessary - please notify (ROE standard)

PFM

22. Metal lingual needed?

Call doctor if metal lingual is needed

23. Metal occlusal needed?

Call doctor if metal occlusal is needed

24. Do you want collars?

*No (ROE standard)

28. Alloy - PFM

NP Selective Laser Melting - Silver

29. Alloy - Full Cast Crown

NP Selective Laser Melting - Silver

Full Denture

30. Denture Stippling

*Yes (ROE standard)

31. Denture Festooning

*Standard (ROE standard)

32. Base and Rims

*Printed denture base (ROE standard)

33. Denture Teeth

*Premium (ROE standard)

34. Set-Ups

Follow Study Model

Splints/Guards

35. Preferred Arch

*Upper arch (ROE standard)

36. Preferred Material

*Hard splint, modeless (ROE standard)

37. Preferred Full Arch Coverage Design

*3.0mm labial coverage, flat plane (ROE standard)

Guided Surgery

38. Always use OEM sleeves?

*No, use non-OEM sleeves when available (ROE standard)

39. Always use sleeveless plastic guides?

*No I like sleeves (ROE standard)

40. Chairside conversion Protocol

*Standard CHROME holes (ROE standard)

41. Always use irrigation ports on plastic guides?

*No (ROE standard)