



CBCT-Based Surgical Guides

800-228-6663 | info@roedentallab.com | www.roedentallab.com
Main Laboratory: 7165 E Pleasant Valley Rd, Independence, OH 44131

Please ✓ applicable boxes. Required fields indicated with asterisk (*). ROE standards (●) apply if no selection is made.

SURGICAL

Arch	Maxillary Mandibular Both
Implant brand	
Implant line	
Fully Guided Kit	
Pilot only	

ROE TO INCLUDE WITH CASE

Fixation Kit (*drill & pins*)
Prosthesis Conversion Luting Material
Abutments
Temporary cylinders
Back-up Denture
Schedule Chairside Assistance

DESIGN


Guide Type	Tooth-supported Tissue-supported Bone-supported CombiGuide (<i>metal reduction guide</i>) Complete guide ● .STL for in-office print Print from Dr's .STL Non-OEM Sleeves ● OEM Sleeves Sleeveless
Options	Without ports ● With ports
Sleeves	GuidedTEMP™ Shell temp Other
Irrigation Ports	
Temp	
Temp Shade	

RECORDS

Impression Type	Digital IOS Polyvinyl / models
IOS Scanner Brand	
CT Scan	Maxillary Patient CT Scan Maxillary Appliance CT Scan Mandibular Patient CT Scan Mandibular Appliance CT Scan

When scanning **dentate** patients, always keep patients open biting on cotton rolls.
If patient wears a denture, ensure denture fits very well and there is NO soft liner.
Place scan markers or 1mm diameter Gutta Percha markers randomly and take two scans: 1. Patient wearing denture in occlusion 2. Denture alone on cotton rolls or packing foam. To learn more, see www.roedentallab.com/recordsFAQ or call.

Please ✓ implant positions below:



Patient Name*	Date	
Placing Doctor*	License #	
Address*		
City*	State*	ZIP*
Email*	Phone*	

DELIVER CASE TO OFFICE

Restoring Doctor*	License #	
Same as placing doctor	Include in online meeting	
Address*		
City*	State*	ZIP*
Email*	Phone*	

Instructions

Signature* License #*

The person signing this work order accepts responsibility for payment and agrees to pay all collection costs including attorney's fees. A 1.5% (18%/yr.) finance charge will be added to all balances due over 30 days.

For the most up-to-date Rx & forms, please visit www.roedentallab.com/forms