

Guide Sleeve Order Form



Misc. Sizes


Qty	Order Number	Int Dia	Height
	Universal		
0	TUBE51-41-35L	4.1	3.5
0	TUBE525-425-4L	4.25	4
0	TUBE57-47-4L	4.7	4
0	TUBE62-5-4L	5	4
0	TUBE6-506-35L	5.075	3.5
0	TUBE625-525-4L	5.25	4
0	TUBE635-535-4L	5.35	4
0	TUBE72-62-35L	6.2	3.5
0	TUBE72-626-35L	6.26	3.5
\$	-		

Doctor Info

	Name
	Address
	City, St. Zip
	Surgery Date:
	Phone
	email address

Order Number Int Dia Height

	Order Number	Int Dia	Height
	Pilot		
0	TUBE22-12-4L	1.2	4
0	TUBE24-14-4L	1.4	4
0	TUBE26-16-4L	1.6	4
0	TUBE27-17-4L	1.7	4
0	TUBE29-19-4L	1.9	4
0	TUBE31-21-4L	2.1	4
0	TUBE34-24-4L	2.4	4
0			



A - Inside Diameter
0.1mm larger than
drill diameter

B - Implant Studio
Options
3mm, 4mm, 5mm

3mm Outside Diameter

Order Number Int Dia Height DRILL

	Order Number	Int Dia	Height	DRILL
	Pilot			
0	TUBE3-2-4L	2	4	1.9mm
0	TUBE3-21-5L	2.1	5	2mm
0	TUBE3-23-5L	2.3	5	2.2mm
0	TUBE3-24-5L	2.4	5	2.3mm
0				

4mm Outside Diameter

Order Number Int Dia Height DRILL

	Order Number	Int Dia	Height	DRILL
	Pilot			
0	TUBE4-21-5L	2.1	5	2mm
0	TUBE4-23-5L	2.3	5	2.2mm
0	TUBE4-24-5L	2.4	5	2.3mm
0	TUBE4-26-5L	2.6	5	2.5mm
0	TUBE4-28-5L	2.8	5	2.7mm
0	TUBE4-29-5L	2.9	5	2.8mm
0				

Universal

0	TUBE4-33-5L	3.3	5	3.2mm
0	TUBE4-35-5L	3.5	5	3.4mm
0	TUBE4-36-5L	3.6	5	3.5mm
0				

5mm Outside Diameter

Order Number Int Dia Height DRILL

Pilot

0	TUBE5-21-5L	2.1	5	2mm
0	TUBE5-23-5L	2.3	5	2.2mm
0	TUBE5-24-5L	2.4	5	2.3mm
0	TUBE5-28-5L	2.8	5	2.7mm
0	TUBE5-29-5L	2.9	5	2.8mm
0				

Universal

0	TUBE5-3-5L	3	5	2.9mm
0	TUBE5-33-5L	3.3	5	3.2mm
0	TUBE5-35-5L	3.5	5	3.4mm
0	TUBE5-39-5L	3.9	5	3.8mm
0				

Guided Systems

Order Number Int Dia System Name Size Implants

0	TUBE65-55-4L	5.5	MIS Guided	ROE MIS Guided	N/A
0	TUBE6-5-35L	5	NEOSS	ROE Nobel Regular (NBR) / Straumann	N/A
0	TUBE5-41-4L	4.1	Nobel	ROE Nobel Narrow (NBN)	3.5mm
0	TUBE6-5-35L	5	Nobel	ROE Nobel Regular (NBR) / Straumann	4.3mm
0	TUBE69-62-4L	6.2	Nobel	ROE Nobel Wide (NBW)	5mm
0	TUBE6-5-35L	5	Straumann	ROE Nobel Regular (NBR) / Straumann	N/A
0	TUBE5-42-45L	4.2	Zimmer	ROE Master A	3.7mm
0	TUBE62-531-45L	5.31	Zimmer	ROE Master B	4.1mm and larger
0	TUBE65-55-4L	5.5	Camlog	ROE Camlog	N/A
0	TUBE527-435-6L	4.35	BioHorizons	ROE BioHorizons Yellow	3.8mm
0	TUBE61-51-6L	5.1	BioHorizons	ROE Biohorizons Green	4.6mm
0	TUBE73-63-6L	6.3	BioHorizons	ROE Biohorizons Blue	5.8mm
0					

Grand Total

Doctor or Company Name:	Signature	Date

Credit Card	
visa / mscd / amex	Exp Date ccv#
# / / /	<input style="width:100%" type="text"/>
Name on Card	
<input style="width:30%" type="text"/>	<input style="width:70%" type="text"/>